

## Exercise 8.3 – Bile Esculin Test


---

### LEARNING OUTCOMES

1. State the principle and procedure of the bile esculin test.
  2. Explain how the bile esculin test is used to distinguish Group D enterococci from other streptococcal species.
- 

The bile esculin test is used in the identification of Group D streptococci, including species of *Enterococcus*. Unlike many Gram-positive bacteria, enterococci grow in 4% bile, which serves as a selective agent in bile esculin medium. This medium also contains esculin, a carbohydrate derivative. While many organisms can hydrolyze or break down esculin, few other than the enterococci are able to do so in the presence of bile. When esculin is hydrolyzed, a molecule called esculin forms and reacts with ferric (iron) ions added to the medium. This reaction causes a black precipitate and darkening of the agar (Figure 8.5).

Recall from earlier experiments that Group D enterococci are gamma hemolytic on blood agar, meaning that they do not act on red blood cells. Therefore, a positive bile esculin reaction is useful for distinguishing these bacteria from other catalase-negative, Gram-positive cocci.



*Figure 8.5: Enterococcus faecalis (left) hydrolyzes esculin in the presence of bile, forming a black complex of esculin and ferric ions while growth of most other Gram-positive cocci is inhibited (right).*

## Exercise 8.3 – Bile Esculin Test

---

### OBJECTIVE

Distinguish Group D *Enterococcus* from other streptococci based on esculin hydrolysis in the presence of bile.

### MATERIALS

- EQUIPMENT: Inoculating loop, incinerator
- MEDIA: Bile esculin agar slants
- CULTURES: *Enterococcus faecalis*, *Staphylococcus aureus*

### PROCEDURE – STUDENTS WORK IN PAIRS

1. Use tape to label each tube with your initials, date, and organism number.
2. Aseptically inoculate bacteria on the slant surface in a single line from bottom to top.
3. Place tubes in a common rack for incubation at 37°C for 18-24 hours.

### FOLLOW UP

1. Evaluate tubes for esculin hydrolysis by observing blackening of the agar.
2. Record results in the lab report.
3. Remove tape from tubes and place them in a common rack for autoclaving.